

INSTRUCTION MANUAL
MODEL 6030/6025
LOW LEVEL AREA MONITOR

Sept 2000
Version 1

Health Physics Instruments
330 South Kellogg Ave, Suite D
Goleta, CA 93117
(805) 967-8422
Division of Far West Technology

Fuse at ½ A Slo Blow 5 x 20 mm

TABLE OF CONTENTS

I.	INTRODUCTION	2
II.	OPERATION	2
III.	INSTALLATION	2
	DETECTOR	2
	DISPLAY	3
	DETECTOR AND DISPLAY CONNECTIONS	3
	I/O CONNECTIONS	3
	CONNECTOR PART NUMBERS	3
	INTERFERENCE	4
	CONTACT CLOSURE FOR EXTERNAL ALARM	4
	RS232	4
	EXTERNAL INTEGRATOR RESET BUTTON	4
IV.	FRONT PANEL	4
V.	ALARMS	5
	ALARM ACTIVATION	5
	ALARM 1	5
	ALARM 2	5
	ALARM 3	5
	FAIL	5
	PRIORITY	6
	NO ALARM	6
	FACTORY SETTINGS	6
	CHECK DISPLAY OR VIEWING THE ALARM CONDITIONS	6
VI.	SETUP	7
	PRESET	8
	ZERO	8
	TRIPSET	8
	DELAY	8
	PAUSE	8
	ALM SETUP (ALARM SETUP)	8
	BUZZER	9
	DECIMAL	9
	START TIME	9
	FAIL TIME	10
	TURNING THE ALARMS OFF	10
VII.	RS232 SERIAL OUTPUT	10
	STATUS BYTE	11
	BIT FUNCTION	11
VIII.	CALIBRATION	12
	ZERO ADJUSTMENT	12
	CALIBRATION PROCEDURE	12
IX.	MAINTENANCE	14
	SHORTCUTS & OPERATIONAL HINTS	14
	ZERO DRIFT	14
	FUSE	14
	LCD CONTRAST	14
	DESICCANT	14
	FAILURE 1	14
	ADJUSTMENTS	15
	MAINTENANCE MODE	16
	TESTING DISPLAY	16
X.	CIRCUIT DESCRIPTION	16
	GENERAL	16
	DETECTOR	17
	DISPLAY	17
XI.	APPENDIX I	18
	USER TABLES	18
XII.	PARTS LIST	20
XIII.	SCHEMATICS	27

I. INTRODUCTION

The Model 6030/6025 Area Monitor is designed to measure low levels of gamma and neutron radiation. It has several alarms and trip points that are user set. The front panel display shows the radiation level and the status of the instrument. The detector is model number 6030 and the display/readout is model number 6025.

II. OPERATION

The operation of the 6030/6025 is very simple. The first line of the display shows the radiation level from the detector. It is updated every second. The second line shows the integrated dose and is updated every second. The lights on the front panel also show the status of the alarms. The front panel controls reset the integrated dose, set the alarm levels, and change the operational parameters.

When the instrument is first turned on there is a wait until the instrument stabilizes. The number of seconds remaining until normal operation begins is shown on the display along with the version number of the software. The detector may take longer than the wait period to completely stabilize to background levels, however this should not effect the alarms because the level is close to background.

III. INSTALLATION

The 6030/6025 is easily installed. It consists of 2 parts, the detector unit (6030) and the display or readout unit (6025). There is a single 4 wire connection between the two and a short cable is supplied with the instrument.

The detector and display can be remoted up to 100 feet. If you wish to use a longer cable you will need to assemble your own cable. The connectors at each end of the cable are crimp type that can be easily assembled in the field. We recommend Belden 8723 cable. This is a 2 pair shielded cable. For short distances any 4 conductor cable will work. For longer distances or if it is used in a noisy environment we recommend a shielded cable. The connectors are wired the same on both ends of the cable.

After setting up the instrument, review the section on ALARMS to know how to change the alarm configurations, on SETUP to know how to change the alarm and zero settings, and on FRONT PANEL to know how to operate the controls on the front panel.

DETECTOR

The detector is installed by first fixing the detector mounting bracket. The detector unit may be mounted in any orientation. It is best if the side of the detector faces the radiation source. The detector sets into the mounting bracket and is held in place with the two thumbscrews. Plug one end of the cable into the circular connector on top of the detector.

DISPLAY

Plug the cable from the detector into the smaller circular connector on the back of the unit. Plug the power cable into the back of the unit and the other end into the wall (110 VAC 60 Hz). Turn the instrument on using the on-off switch on the back panel of the display. The display should turn on and show the warm-up time counting down.

DETECTOR AND DISPLAY CONNECTIONS

Small circular connector on the back of the display and on the top of the detector.

PIN #	DESCRIPTION	WIRE COLOR To Display
1	DATA	WHITE
2	DATA	GREEN
8	+12 VDC	RED
9	GROUND	BLACK AND SHIELD

I/O CONNECTIONS

Large circular connector on the back of the display.

PIN #	DESCRIPTION
1	RELAY N.C. (opens on alarm and power off)
2	RELAY COMMON
3	RELAY N.O. (closes on alarm and power off)
4	EXTERNAL RESET SWITCH
8	TXD FOR RS232
9	RXD FOR RS232 (Not used)
10	GROUND FOR RS232
13	+12 VOLTS (not used)
14	GROUND for External Reset Switch

CONNECTOR PART NUMBERS

This table lists the connectors that mate with the connectors in the 6030/6025. All are plastic and made by AMP.

PART	SMALL CIRCULAR CONNECTOR	LARGE CIRCULAR CONNECTOR
CONNECTOR PLUG	AMP 20678081 13-9 std	AMP 206044 17-14 rev
CABLE CLAMP	AMP 206966-1	AMP 206070-1
CONTACTS	AMP 661094 (socket)	661074 (pin)

INTERFERENCE

Interference can cause annoying false alarms. The most likely cause of interference is from noisy AC switches and lamp dimmers. Plugging the offending product into a line filter can often eliminate the interference.

CONTACT CLOSURE FOR EXTERNAL ALARM

The contact closure from the internal relay is on pins 1,2 and 3 of the large circular connector on the rear of the display. Normally an external alarm would be connected between pin 2 and 3. This alarm operates in the fail-safe mode. If you were to remove power to the display, the relay would open closing contacts on pin 2 and 3. If you are using this relay in a interlock system, connect it between pin 1 and 2. Contacts rated 0.5 A 115 VAC, 1 A 24 VDC resistive.

RS232

An RS232 output is available between pin 9 (TXD) and pin 10 (ground) of the large circular connector on the rear of the display. The RXD (pin 8) is not supported with this software.

EXTERNAL INTEGRATOR RESET BUTTON

If you want an external reset button for the integrate range, connect a N.O. switch between pin 4 and 14 on the large circular connector on the back of the display. Push the button when you want to reset the integrated dose. You should hold the button down for 1 second to assure that the integrate range is reset. To test this connection see Maintenance Mode.

IV. FRONT PANEL

The front panel of the display contains the LCD display, the alarm LEDs, the beeper, and the four pushbuttons.

The LCD shows the radiation level on the top line and the integrated dose on the bottom line. There is a dot on the right hand side of the display that flashes every second to show that the instrument is functioning. There is no overrange indication. The maximum level is shown on the calibration report.

There are 4 LED lights along the bottom of the front panel. Their action and colors are from left: unused, red-alarm1, yellow-alarm2, and green-no alarms and no fail. The factory settings use the red LED for the alarm, the yellow LED for a warning and the green LED to indicate that all is OK.

The 4 pushbuttons control the setting of the alarms, the fail and turn on delays, and resetting the integrated dose. Pushing the MENU button will cycle through the two possible settings. The first display allows resetting the integrated dose and the second allows setting the alarm levels and other alarm and operational parameters.

To reset the dose, push MENU until the display shows "Integrated Dose?" Then push both YES buttons at the same time. The settings menu works the same way. Push MENU until the display shows "SETTINGS?" Then push both YES buttons at the same time. See the section on ALARMS and SETTINGS for more information.

V. ALARMS

There are 4 alarms built into the instrument. All are based on the radiation level and all are adjustable both for the level and for their control over the indicators and relays. The only thing that is not adjustable is their priority.

Alarm 1, 2, and 3 all trip if the count level exceeds the alarm setting. The Fail alarm occurs if there are no counts from the detector for a preset time. Each alarm has different settings, and each alarm is designed to look at a different part of the level. Each alarm may be individually turned off if it is not needed by setting the trip level to 0. Alarm 1, 2, and 3 all have several settings associated with them. The different parts of the settings are:

1. **TRIP SET** This is the alarm level. This setting is a number that is compared with the level from the detector.
2. **DELAY** This is the number of intervals that the alarm must be consecutively activated to actually trip the alarms.
3. **PAUSE** This is the number of seconds after the level has decreased below the trip set that the alarm will remain activated. It is used to keep the alarm on longer than one second. It is usually set around 10 seconds.

ALARM ACTIVATION

Alarms 1, 2, and 3 work in the following manner. When the level rises above the tripset, the delay counts down every interval period from its preprogrammed level. When it reaches zero it turns the alarm on. Until the alarms are activated if the level decreases below the tripset, the delay will reset to its preprogrammed level. This helps to keep noise from tripping the alarm. The pause works like the delay in that it too counts down every second from its preprogrammed level once the level drops below the tripset. When it reaches zero the alarms will be deactivated for that alarm. If other alarms are still activated they will continue to be activate. If during the pause period the level rises above the tripset again, even for one interval, the pause will be reset to the preprogrammed level. Thus once activated the alarms will stay on for at least the pause period following the last occurrence of a trip. This helps to keep the alarms from cycling on and off in a marginal situation.

ALARM 1

This alarm has the highest priority. It trips if the level rises above the tripset. It is usually setup to trip the front panel RED LED and to activate the relay and beeper. When tripped it will show ALARM 1 on the second line of the CHECK DISPLAY. It is usually set with the highest alarm level, a short delay and a 10 second pause.

ALARM 2

This alarm is identical to alarm 1 except it has the second highest priority. It is usually setup to trip the front panel YELLOW LED and to not activate the relay and beeper. It is used with a short delay and a 10 second pause. When tripped it will show ALARM 2 on the second line of the check display.

ALARM 3

This alarm has the third highest priority. It is usually not activated but can be used for any additional configuration. For example, it could be set to a very low level with a long delay to look for slowly changing levels. When tripped it will show ALARM 3 on the second line of the CHECK DISPLAY.

FAIL

This alarm is used to indicate that the detector is not functioning. It turns off the alarms and indicators and shows FAILURE in the display. It will not activate the relay because there is no known hazard. This alarm will activate if there are no counts from the detector for a number of seconds. The number of seconds is the only setting. When tripped it will show FAIL on the second line of the CHECK DISPLAY.

PRIORITY

The following table shows the priority of the alarms. Alarm 1 has the highest priority which means it will supersede the lower priority alarms. Fail has the lowest priority because if the instrument fails, it obviously has no counts and cannot set the other alarms.

- Alarm 1 Highest Priority
- Alarm 2
- Alarm 3
- Fail Lowest Priority

NO ALARM

The NO ALARM condition is usually used with the GREEN LED turned on, the relay and beeper turned off and the display showing OK on the second line of the CHECK DISPLAY.

FACTORY SETTINGS

The following table lists the factory settings for the alarms.

ALARM	USE	DELAY	PAUSE	RED LED	YELLOW LED	GREEN LED	RELAY	BEEPER
1	High Levels	0	10	ON	OFF	OFF	ON	ON
2	Warning	0	10	OFF	ON	OFF	OFF	OFF
3	Not used							
FAIL	Detector Failure	N/a	N/a	OFF	OFF	OFF	OFF	OFF
Normal	No Trip	N/a	N/a	OFF	OFF	ON	OFF	OFF

CHECK DISPLAY OR VIEWING THE ALARM CONDITIONS

The CHECK DISPLAY shows the pause and trip status on the 2nd line of the LCD during normal operation by pushing the right hand push-button. When the push-button is down, the status of the three alarms will be displayed on the LCD.

When both right hand buttons are pushed at the same time, then the display changes to three groups of characters/digits. Each represents an alarm and are in the following order from left to right: alarm 1, alarm 2, alarm 3 and the fail time. If an alarm is off its values will be displayed as *--.

The first character of each group is an '*' if the alarm is not tripped and a 'T' if the alarm is tripped. The next number is a 2 digit hex number of either the delay or the pause. If the alarm is not tripped then it is the delay. If the alarm is tripped then it is the pause.

This example will use alarm 1. Assume that the delay is set to 5 and the pause is set to 8. Normally, with no trip it would read '*05'. The '*' indicates it is not tripped and the '05' is the delay. If the level was brought higher than the tripset then the delay would start to count down every interval period until it reached zero. This shows the delay period. If the level were to decrease below the trip set during the time it was counting down, then the delay would revert back to its level which is 5. When the delay reaches zero, the indicator will change from a '*' to a 'T' to indicate it has been tripped, and the alarms will be set. The display will then show 'T8' and will continue to show 'T8' until the level is brought down below the tripset. When the level is brought below the tripset, the pause will start counting down, decreasing by 1 every second. When it reaches zero, the 'T' will change back to a '*' and the alarms will be set to the no alarm condition.

The fail time is a hex value of the seconds remaining before the fail mode will be activated. It will start counting down if there is a failure. When it reaches zero it sets the alarms to the fail setting.

VI. SETUP

The instrument has been setup with its preset values. These values are programmed into the EEPROM (changeable permanent memory). They can be changed by the user. This section shows how to change these presets. APPENDIX I contains blank forms for recording your settings. We recommend that you copy this page and use it to figure out your changes.

The setup mode is different from the normal operation of the instrument. To enter into the setup mode first push MENU twice then the two YES buttons at the same time. **NOTE:** To exit this mode, keep pushing the MENU button until you have cycled through all of the choices.

The buttons on the control panel will do the following:

MENU (left hand button) will bring up the next item to adjust. Repeatedly pushing the MENU button will cycle through all the adjustments.

RESET (left of center button) is not used.

UP(right of center button) will move the arrow on the bottom line from one digit to the next. Every time it is pushed the arrow will move left to the next digit. When it gets to the last digit it will jump to the first digit.

DOWN (far right hand button) will increment the digit that the arrow points to. Every time the button is pushed the digit will increase.

It only takes a few seconds of playing with the buttons to understand how they function.

Some adjustments have 3 digits and some have 4. The adjustments with 3 digits have a maximum setting of 255. If they are set above 255 they will actually be set to 255. The bottom line of the display reminds you that they have a maximum value of 255. The four digit adjustments have no restrictions, they can be adjusted from 00.00 to 99.99.

Repeatedly pushing the MENU button will cycle the display through all of the adjustments. After the last adjustment the program will go back to normal operation. Most of the settings are saved in EEPROM after the last item which is the fail-safe time, consequently if you are part way through changing the settings and decide you don't want the new values you can turn power off then back on or if the top cover is removed push the reset button (S1) which is under the ribbon cable to the display.

The following is a list of the parameters in the order that they are seen on the display along with the factory presets. The letters A1, A2 etc. refer to alarm 1, alarm 2 etc. DELAY A2 is the delay value for alarm 2.

Presets @255	
ZERO	60
TRIPSET 1	1.00
TRIPSET A2	0.80
TRIPSET A3	0
DELAY A1	1
DELAY A2	1
DELAY A3	10
PAUSE A1	10
PAUSE A2	10
PAUSE A3	10
Alm Setup A1	137
Alm Setup A2	137
Alm Setup A3	4
No alm Setup	2
Start Time	255
Fail Time	100

Please read the section on the alarms to become familiar with the action of the alarms. All of the parameters are reviewed below. You must cycle through all of the parameters to get back to a normal display.

PRESET

Set this value to 255 if you want all of the adjustments to be reset to their factory preset values. If you do not want the factory preset settings, then push MENU again to go to the next item.

ZERO

This is the value that is subtracted from the reading from the detector. The level in the detector is elevated above zero. See the calibration section.

TRIPSET

This is the alarm level. This setting is a 4 digit number that is compared with the reading from the detector.

DELAY

This is the number of intervals that the alarm must be consecutively activated to actually trip the alarms.

PAUSE

This is the number of seconds after the level has decreased below the trip set that the alarm will remain activated. It is used to keep the alarm on longer than one interval. It is usually set around 10 seconds. It can be set longer but it usually is determined by how long the beeper (or external alarm) needs to be on to arouse someone that there is a problem.

ALM SETUP (ALARM SETUP)

This is a number that is used to set the condition of the alarms, indicators and external outputs. There are 4 setups, one for each of the three alarms and one for no alarms. Below is a description of the alarms, indicators and external outputs. Each can be set to only two values, 0 or 1. Following the descriptions is the method used to calculate the values and to determine the decimal value.

RELAY

This controls the relay. The contacts of the relay are brought out to the large circular connector. When the relay is set to 0, the relay is closed. When the relay is 1, the relay is open. The relay is operated in the fail-safe mode (OFF actually energizes the relay). The signal to control the relay is also routed to pin 8 on P4 on the main circuit board in the display.

P4:5

This is pin 5 on plug P4 which is located on the main circuit board in the display. When the relay is set to 1 then this pin is high. This pin has no normal function and is not used in normal operation of the monitor.

P4:6

This is pin 6 on plug P4 which is located on the main circuit board in the display. When P4:6 is set to 1 then this pin is high. This pin has no normal function and is not used in normal operation of the monitor.

P4:7

This is pin 7 on plug P4 which is located on the main circuit board in the display. When P4:7 is set to 1 then this pin is high. This pin has no normal function and is not used in normal operation of the monitor.

P4:8

This is pin 8 on plug P4 which is located on the main circuit board in the display. When P4:7 is set to 1 then this pin is high. This pin has no normal function and is not used in normal operation of the monitor.

RED LED

This is the front panel red LED. If the RED LED is set to 1 then the LED is on.

YELLOW LED

This is the front panel yellow LED. If the YELLOW LED is set to 1 then the LED is on.

GREEN LED

This is the front panel green LED. If the GREEN LED is set to 1 then the LED is on.

BUZZER

This is the front panel buzzer. If BUZZER is set to 1 then the buzzer is turned on and emits a loud continuous beep.

DECIMAL

The following table is a compilation of the settings of all the parts of the alarm setup. The value is a decimal value that is calculated from the results of the table. The line of one's and zero's on a row is actually a binary number. This number is converted to decimal and that is the decimal number. For example the first line of the table below shows:

alarm * * * 0 11 0 0

If you change the * to zeros it becomes the number 00001100. This is a binary number. To find its decimal equivalent, look at the binary to decimal conversion table in the appendix. Look at the first binary column from the left and about 13 numbers down. You should find the number 00001100. Next to it is the number 12. This is the decimal conversion. 00001100 in binary is 12 in decimal. This decimal number is the number you enter into the alarm setup.

The following is a table that shows the normal operation of the instrument as it is setup using the factory presets.

	P4:8	P4:7	P4:6	RELAY and P4:5	BEEPER	RED LED	YELLOW LED	GREEN LET	DECIMAL VALUE
Alarm 1	*	*	*	0	1	1	0	0	12
Alarm 2	*	*	*	1	0	0	1	0	18
Alarm 3	*	*	*	1	0	0	0	0	16
No Alarm	*	*	*	1	0	0	0	1	17
FAIL	0	0	0	1	0	0	0	0	16

- *Setting does not matter for normal operation. We suggest each of these be set to 0. The decimal calculations assume that the items marked * are set to 0. FAIL is only show for information purposes, it is not changeable.

START TIME

When the instrument is turned on it waits before going into normal operation. This gives time for the instrument to stabilize. The display shows the seconds counting down until normal operation. The start time is the starting number for the countdown. It can also be thought of as the start delay. If it is set to a low number the detector may not stabilize in time causing the alarms to trip.

FAIL TIME

If the display does not receive a signal from the detector it will show a failure in the display. The fail time is the number of seconds after the last signal before the instrument will activate the failure mode. If the fail time is set to 60 seconds, it will take 60 seconds after the detector fails before the instrument will display FAILURE. The factory setting is 30. The detector should send a signal representing 0.01 mR/h at least every second if it is operating normally.

TURNING THE ALARMS OFF

Each of the alarms can be turned off. To turn off alarm 1, 2 or 3, set the tripset to zero.

VII. RS232 SERIAL OUTPUT

The serial output only sends RS232 data. It does not receive data.

SIGNON MESSAGE

At turn-on the instrument sends a signon message then a packet of data. The message is:

```
HPI 6030 VER 1.0  
000064 000050 000000 00 00 00 00 00 00
```

The version number of the software is shown on the first line.

The data on the second line from left to right:

1	XXXXXX	3 hex digits Alarm trip level for alarm 1
space		
2	XXXXXX	3 hex digits Alarm trip level for alarm 2
space		
3	XXXXXX	3 hex digits Alarm trip level for alarm 3
space		
4	XX XX	2 hex digits separated by a space. The first is the delay. The second is the pause. Both are for alarm 1.
space		
5	XX XX	2 hex digits separated by a space. The first is the delay. The second is the pause. Both are for alarm 2
space		
6	XX XX	2 hex digits separated by a space. The first is the delay. The second is the pause. Both are for alarm 3

SERIAL STREAM

Every second the display sends out a packet of data via the RS232 serial output. The packet is the data at the time it was sent and should look like:

```
000001 0000000000 00
```

the data from left to right is:

1	XXXXXX	3 hex digits Counts per second from detector. The same value as on the display. When converted to decimal add decimal point 2 digits from the right: X.XX mR/h.
---	--------	---

- 2 XXXXXXXX 5 hex digits integrated dose. When converted to decimal add decimal point 2 digits from the right: X.XX mR.
- 3 XX 1 hex digit This is the status byte. See below for a description.

STATUS BYTE

The status byte consists of 8 bits. The bits are represented as follows: Bit 0 is the LSB and bit 7 is the MSB. Only bits 0 thru 4 are used. Use the lookup table in the appendix to change from hex to binary.

BIT FUNCTION

- 0 0=no trip, 1=trip for alarm 1
- 1 0=no trip, 1=trip for alarm 2
- 2 0=no trip, 1=trip for alarm 3
- 3 0=no fail, 1=fail

VIII. CALIBRATION

The instrument has three adjustments for calibration, two are for the zero and the other is the calibration adjust.

ZERO ADJUSTMENT

The detector is biased above zero, i.e. when it is exposed to no radiation it still sends a signal to the display. The display digitally subtracts a zero offset value from this level. This is the value that is shown on the LCD. It is possible to show the level without the zero offset subtracted from the reading by pushing the right hand button and holding it down. The level on the top line is the level without the zero subtracted. The right hand number on the bottom line is the zero offset value. It can be changed through the setup menu. The factory setting is 0.60 mR/h.

There are two ways to approach the zero adjustment. The first is to adjust the zero offset in the SETUP menu to the factory setting of 0.60 mR/h then adjust the zero trimmer in the detector until the normal display shows a value above zero, then to slowly back it off to zero. The second is to leave the zero adjustment in the detector where it is and adjust the zero offset in the display until the display shows zero. There is some noise on the signal and either adjustment is made to just eliminate any reading when the detector is not exposed to any radiation. It is a good idea to monitor the zero for a few minutes to make sure it is correct.

We recommend a combination of the two. First set the zero offset in the display to 0.60 mR/h. Then adjust the zero trimmer in the detector until the display just reads zero most of the time. Then use the zero offset as a fine adjustment if you need to. The zero offset in the display should be set somewhere from 0.30 to 0.80 mR/h. If it is set too low then the fail may not work properly. If it is set too high then the reading may drift. Remember that if the instrument is reset to the factory settings then the zero offset will be set to 0.60 mR/h.


You can check the zero adjustment by integrating for several hours in the Maintenance Mode. The integrated dose (mR) will show you where the zero is set because Maintenance Mode removes the zero offset from both the rate and integrate ranges. It is always best to let the instrument sit for several hours before adjusting the zero.

CALIBRATION PROCEDURE

Figure 1 shows the locations of the adjustments on the circuit board in the detector housing.

1. Connect the display to the detector and turn the instrument on. Wait for the stabilization period to end.
2. Check and adjust zero by adjusting either the zero in the detector or the zero offset in the display or both. See the ZERO ADJUSTMENT section above.
3. Expose the detector to a known quantity of radiation. Adjust the calibration adjust in the detector until the reading in the CHECK DISPLAY is the correct value.
4. Check zero and readjust if necessary. Alternate between step 2 and 3 until both zero and the calibration are correct.
5. Check at least 1 point on every decade. (1, 10, 100 mR/h)
6. Continue to expose the detector to higher and higher levels until the display reaches a maximum, usually from 100 to 300 mR/h. Note this maximum level on the calibration report.
7. The integrate range should be checked. Do not exceed the maximum dose rate. There are 2 ways to check this range. The first is to turn time the source and the second is to time the exposure after the instrument is reading the correct rate. For the first method of timing the source, wait for the exposure rate to drop to zero before starting the exposure (and timer). After the timer has timed out and exposure stopped, wait for the exposure rate on the display to return to 0 before taking the exposure numbers from the display. You want the ramp up and down.

For the second method, expose the instrument to the calibrated rate and wait for the rate display to stabilize. Reset the integrated dose and start your timer. When the time is up, read the dose on the fly.



DETECTOR

Figure 1 Calibration And Zero Adjustment Locations

IX. MAINTENANCE

This section discusses the circuit of the instrument and any adjustments that may be needed. It also describes a test mode for checking out the instrument.

SHORTCUTS & OPERATIONAL HINTS

The 6030/25 has several shortcuts.

Push the RESET button on the front panel of the display during the warm-up period to cancel the warm-up period.

Pushing the right hand button during normal operation to see the status and to see the zero offset.

Pushing the two right hand buttons together during normal operation to show the status of the alarms and fail.

Push the RESET button on the control panel during normal operation to cancel the local buzzer.

Push the DOWN BUTTON and keeping it on during power up will enter the maintenance mode.

If you have the complete cover off of the display, the reset button (S1) on the top of the circuit board can be used to abort the setup routine. Just push it while in the setup routine. You can also abort the setup routine by turning off the power.

If you are testing the alarms and the noise is too loud, put a piece of tape over the beeper. It will not make it quiet but it will reduce the volume.

ZERO DRIFT

If you notice the zero drifting change the zero offset in the SETUP Menu to compensate. Remember you want the reading on the display to read zero most of the time if there is only background radiation. You can see where the zero is by looking at the CHECK DISPLAY. It is best to let the instrument warm up for a day or so before adjusting zero.

FUSE

The fuse is located in the socket for the line cord. Remove the line cord and pry out the fuse holder. Replace the fuse with a 5 x 20 mm fuse rated at ½ A slow blow. Fast-acting will also work in an emergency however we have found that they occasionally blow because of the type of load.

LCD CONTRAST

The LCD contrast is on the large circuit board in the display. Remove the top cover and turn R2 to adjust the contrast. It is located under the ribbon cable that goes to the LCD display .

DESICCANT

The desiccant is inside the electrometer housing in the detector box. It only needs to be replaced if the zero is excessively high or erratic.

FAILURE 1

If the display should show FAILURE 1 this indicates that the eeprom is defective and that all of the settings including the alarm levels and zero offset may be corrupted. You can try to re-program the eeprom to the factory settings by pushing MENU. This will change the display to the setup menu. Then try setting the preset settings. If it still shows FAILURE 1 then replace the eeprom and repeat with the factory settings.

ADJUSTMENTS

The zero and gain adjustments for the detector are covered in the calibration section.

HIGH VOLTAGE ON DETECTOR CIRCUIT BOARD


The high voltage for the detector is adjusted on the circuit board with R37 to -800 volts. Use a high impedance voltmeter ($>10E12$ ohm) to measure the voltage at the junction of C27 and R30.

NEGATIVE VOLTAGE ON DETECTOR CIRCUIT BOARD

The negative voltage for the preamp electrometer is adjusted to 0.1 volt. Place a voltmeter between ground and pin 6 of U9. Adjust R26 until the voltage is 0.1 volt.

POWER SUPPLY READJUSTMENT

The 12 VDC 110/220 Volt power supply located under the circuit board in the display is readjusted for an output of 11 volts. This reduces the power dissipation of the regulators on the display and detector boards.



DETECTOR



DISPLAY/READOUT

Figure 2 Adjustment Locations

MAINTENANCE MODE

This is a special display that allows you to monitor the rate, integrate, maximum and minimum level of the detector.

Use this mode to check the zero reading from the display. You can enter this mode and let it sit for a day or two to check both the maximum and minimum level. This shows if there is excessive noise in the display. The integrated dose can also show if the average is too high. Remember to mentally subtract the zero offset from the readings, both rate and integrate.

To enter into Maintenance Mode, hold down the DOWN button while you turn on the display.

The display shows the radiation level on the top line of the LCD followed by the maximum level since reset. The bottom line shows the integrated dose followed by the minimum level since reset. The radiation level, and minimum and maximum levels all read in mR/h. The integrated dose is in mR.

All levels are actual levels without the zero offset subtracted from the readings.

When the external integrate reset button is pushed, the lower right hand digit will change to a * .

Pushing RESET will reset the minimum and maximum levels.

Pushing the UP button will activate the relay.

Pushing the MENU button will cycle through the alarms and annunciators.

TESTING DISPLAY

The display alarms and annunciators can be tested by pushing down the MENU button on the front panel. When it is held down the instrument will cycle through 10 different annunciators in the following

order:

- Green LED
- Yellow LED
- Red LED
- Front panel buzzer
- Relay and P4:5
- All off
- P4:6
- P4:7
- P4:8
- All lights on

The cycle will then repeat as long as the MENU button is held down. (P4 is not available except on the circuit board.)

X. CIRCUIT DESCRIPTION

GENERAL

The detector sends a frequency to the display that depends on the radiation level. The frequency can vary from 0 to about 20,000 Hz. Each Hz represents 0.01 mR/h. The zero of the electrometer is also biased up from absolute zero. Thus with no radiation the detector has a signal of 0.60 mR/h. A frequency of 360 Hz from the detector is thus displayed as 3.00 mR/h; 360 is 3.60 mR/h and subtract the zero offset of .60 results in 3.00. The signal between the detector and the display is based on a RS485 driver and is a complimentary signal, i.e. both signal lines change state.

DETECTOR

The detector measures a small amount of current from the ion chamber and creates a frequency output. The detector needs a high voltage supply to function properly.

The detector V1 is a 3 liter three terminal ion chamber. Ionizing radiation produces ion pairs inside the detector that are collected on the center signal portion of the detector. This current is integrated onto U101, an electrometer amplifier with a high value feedback resistor (R101). The negative ion chamber current makes the output of the electrometer move positive. The input filter C104 and R 103 smooth any peak transients from pulsed radiation.

The output of the electrometer goes from 0 volt to 4 volts. U1 acts as a buffer amplifier and it buffers the voltage to frequency converter U2. The output of U2 goes to U4:A which is a pulse generator. The pulse from U4:A is sent to the 75176 line driver and then to the cable to the readout.

The high voltage power supply is a switching type generating -800 volts. U10 is an oscillator that feeds the one-shot U10B. The output of this one-shot turns on Q2 creating a pulse in the primary of T1. The secondary of T2 feeds a voltage quadrupler. The high voltage is measured by U11:A which turns the one-shot on and off to regulate the voltage. R48, R30 and C27 and C37 form a filter to remove any ripple from the high voltage. U11:B monitors the high voltage and turns off the output signal if the HV falls too low.

DISPLAY

The display measures the frequency from the detector, subtracts background using software, and sets the alarms and relays according to the settings in memory. It contains the line power supply to power the detector and itself.

The display is controlled by the microprocessor, U3. U1 is an address latch that separates the address and data. U5 is the LCD display. U2 is the EEPROM that stores the variables. U8 is the address decoder for the LCD and for the input latch U19 and the output latch U18. The relay is connected directly to the microprocessor through U9:D and Q1. This is done to make the relay fail-safe. If the microprocessor is reset either by power on or by the watchdog timer U4, then the output pin P33 (~INT) will float high opening the relay. U7 is a space for extra RAM should it be needed.

U12 and U13 are RS232 and RS485 outputs respectively.

U18 is an output latch that feeds the output buffers U16 and U17 this drives the front panel LED's. U19 is an input latch for the pushbuttons and for the external reset switch. U21 is an A/D converter. U24 is a DAC.

The incoming pulses from the detector go to U26 and then to the microprocessor U3 where they are counted.

R1 adjusts the contrast. U4 is a watchdog timer for the MPU and it needs a pulse on P1.0 (clock) at least every second to keep the watchdog from timing out and resetting the microprocessor.

U10 is a 5 volt regulator and U11 generates a negative 5 volts for the LCD.

U20 is an optocoupler that separates the external reset signal from the internal circuit board. It is used to reset the integrator.

MODIFICATIONS

Not all of the parts on the circuit board and schematic are included in this model. They are included for future changes or to allow this circuit board to be used in different ways. The detector board has a place for U14, a microprocessor. The designators marked A on the circuit boards are the same as designators marked without an A.

The display has more options. It includes jumpers for RS232 and RS485 inputs and outputs on both the communications lines and the data lines from the detectors. The lines from the detector are designed to accept pulses or serial data from the detector. The watchdog can be disabled by removing JP1. This board also contains an ACD, a DAC, and extra outputs and inputs that are not used. You may want to consult the factory if you need any of these enabled or changed because several are software dependent.

XI. APPENDIX I

USER TABLES

The blank tables below are useful to determine how you want the alarms setup. Copy it and fill it out not only for ease of programming, but also for your records.

DATE _____ SERIAL NUMBER _____

LOCATION _____ BY _____

ALARM SETUP

	P4:7	P4:6	P4:5	RELAY	BEEPER	RED LED	YELLOW LED	GREEN LET	DECIMAL VALUE
Alarm 1									1
Alarm 2									2
Alarm 3									3
No Alarm									4

ENTER THE DECIMAL NUMBERS FROM THE DECIMAL VALUE BOXES ABOVE INTO THE SETUP BOXES BELOW

SUMMARY OF SETTINGS

	TRIPSET	DELAY	PAUSE	SETUP
ALARM 1				1
ALARM 2				2
ALARM 3				3
NO ALARM				4

START TIME _____

FAIL TIME _____

.APPENDIX II

DECIMAL/HEX/BINARY
 CONVERSION TABLE

DEC HEX BINARY

0 00 00000000
 1 01 00000001
 2 02 00000010
 3 03 00000011
 4 04 00000100
 5 05 00000101
 6 06 00000110
 7 07 00000111
 8 08 00001000
 9 09 00001001
 10 0A 00001010
 11 0B 00001011
 12 0C 00001100
 13 0D 00001101
 14 0E 00001110
 15 0F 00001111
 16 10 00010000
 17 11 00010001
 18 12 00010010
 19 13 00010011
 20 14 00010100
 21 15 00010101
 22 16 00010110
 23 17 00010111
 24 18 00011000
 25 19 00011001
 26 1A 00011010
 27 1B 00011011
 28 1C 00011100
 29 1D 00011101
 30 1E 00011110
 31 1F 00011111
 32 20 00100000
 33 21 00100001
 34 22 00100010
 35 23 00100011
 36 24 00100100
 37 25 00100101
 38 26 00100110
 39 27 00100111
 40 28 00101000
 41 29 00101001
 42 2A 00101010
 43 2B 00101011
 44 2C 00101100
 45 2D 00101101
 46 2E 00101110
 47 2F 00101111
 48 30 00110000
 49 31 00110001
 50 32 00110010
 51 33 00110011
 52 34 00110100
 53 35 00110101
 54 36 00110110
 55 37 00110111
 56 38 00111000
 57 39 00111001
 58 3A 00111010
 59 3B 00111011
 60 3C 00111100
 61 3D 00111101
 62 3E 00111110
 63 3F 00111111

64 40 01000000
 65 41 01000001
 66 42 01000010
 67 43 01000011
 68 44 01000100
 69 45 01000101
 70 46 01000110
 71 47 01000111
 72 48 01001000
 73 49 01001001
 74 4A 01001010
 75 4B 01001011
 76 4C 01001100
 77 4D 01001101
 78 4E 01001110
 79 4F 01001111
 80 50 01010000
 81 51 01010001
 82 52 01010010
 83 53 01010011
 84 54 01010100
 85 55 01010101
 86 56 01010110
 87 57 01010111
 88 58 01011000
 89 59 01011001
 90 5A 01011010
 91 5B 01011011
 92 5C 01011100
 93 5D 01011101
 94 5E 01011110
 95 5F 01011111
 96 60 01100000
 97 61 01100001
 98 62 01100010
 99 63 01100011
 100 64 01100100
 101 65 01100101
 102 66 01100110
 103 67 01100111
 104 68 01101000
 105 69 01101001
 106 6A 01101010
 107 6B 01101011
 108 6C 01101100
 109 6D 01101101
 110 6E 01101110
 111 6F 01101111
 112 70 01110000
 113 71 01110001
 114 72 01110010
 115 73 01110011
 116 74 01110100
 117 75 01110101
 118 76 01110110
 119 77 01110111
 120 78 01111000
 121 79 01111001
 122 7A 01111010
 123 7B 01111011
 124 7C 01111100
 125 7D 01111101
 126 7E 01111110
 127 7F 01111111

128 80 10000000
 129 81 10000001
 130 82 10000010
 131 83 10000011
 132 84 10000100
 133 85 10000101
 134 86 10000110
 135 87 10000111
 136 88 10001000
 137 89 10001001
 138 8A 10001010
 139 8B 10001011
 140 8C 10001100
 141 8D 10001101
 142 8E 10001110
 143 8F 10001111
 144 90 10010000
 145 91 10010001
 146 92 10010010
 147 93 10010011
 148 94 10010100
 149 95 10010101
 150 96 10010110
 151 97 10010111
 152 98 10011000
 153 99 10011001
 154 9A 10011010
 155 9B 10011011
 156 9C 10011100
 157 9D 10011101
 158 9E 10011110
 159 9F 10011111
 160 A0 10100000
 161 A1 10100001
 162 A2 10100010
 163 A3 10100011
 164 A4 10100100
 165 A5 10100101
 166 A6 10100110
 167 A7 10100111
 168 A8 10101000
 169 A9 10101001
 170 AA 10101010
 171 AB 10101011
 172 AC 10101100
 173 AD 10101101
 174 AE 10101110
 175 AF 10101111
 176 B0 10110000
 177 B1 10110001
 178 B2 10110010
 179 B3 10110011
 180 B4 10110100
 181 B5 10110101
 182 B6 10110110
 183 B7 10110111
 184 B8 10111000
 185 B9 10111001
 186 BA 10111010
 187 BB 10111011
 188 BC 10111100
 189 BD 10111101
 190 BE 10111110
 191 BF 10111111

192 C0 11000000
 193 C1 11000001
 194 C2 11000010
 195 C3 11000011
 196 C4 11000100
 197 C5 11000101
 198 C6 11000110
 199 C7 11000111
 200 C8 11001000
 201 C9 11001001
 202 CA 11001010
 203 CB 11001011
 204 CC 11001100
 205 CD 11001101
 206 CE 11001110
 207 CF 11001111
 208 D0 11010000
 209 D1 11010001
 210 DA 11010010
 211 D3 11010011
 212 D4 11010100
 213 D5 11010101
 214 DE 11010110
 215 D7 11010111
 216 D8 11011000
 217 D9 11011001
 218 DA 11011010
 219 DB 11011011
 220 DC 11011100
 221 DD 11011101
 222 DE 11011110
 223 DF 11011111
 224 E0 11100000
 225 E1 11100001
 226 E2 11100010
 227 E3 11100011
 228 E4 11100100
 229 E5 11100101
 230 E6 11100110
 231 E7 11100111
 232 E8 11101000
 233 E9 11101001
 234 EA 11101010
 235 EB 11101011
 236 EC 11101100
 237 ED 11101101
 238 EE 11101110
 239 EF 11101111
 240 F0 11110000
 241 F1 11110001
 242 F2 11110010
 243 F3 11110011
 244 F4 11110100
 245 F5 11110101
 246 F6 11110110
 247 F7 11110111
 248 F8 11111000
 249 F9 11111001
 250 FA 11111010
 251 FB 11111011
 252 FC 11111100
 253 FD 11111101
 254 FE 11111110
 255 FF 11111111

XII. Parts List

The designators marked A on the circuit boards are the same as designators marked without an A.

DISPLAY READOUT

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
R1	1		10K	Resistor Network, 1/4W, 5%, CF		SBC3-001 Display
R2	1		20K	Trimmer, 1T		SBC3-001 Display
R3	1		10K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R4	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R5	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R6	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R7	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R8	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R9	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R10	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R11	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R12	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R13	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R14	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R15	0		0 Ohm	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R16	1		0 Ohm	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R17	1		0 Ohm	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R18	1		0 Ohm	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R19	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R20	1		10K	Resistor network, 1/4W, 5%, CF		SBC3-001 Display
R21	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R22	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R23	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R24	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R25	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R26	1		1K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R27	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R28	1		10K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R29	1		510	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R30	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R31	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R32	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R33	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R34	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R35	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R36	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R37	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R38	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R39	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R40	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R41	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R42	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R43	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R44	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R45	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R46	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R47	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R48	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
R49	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R50	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R51	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R52	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R53	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
R54	1		100K	Resistor, 1/4W, 5%, CF		SBC3-001 Display
R55	0			Resistor, 1/4W, 5%, CF		SBC3-001 Display
C1	1		22 pF	Capacitor, mono ceram		SBC3-001 Display
C2	1		22 pF	Capacitor, mono ceram		SBC3-001 Display
C3	1		0.47 uF	Capacitor, mono ceram		SBC3-001 Display
C4	1		470 uF	Capacitor, Electro		SBC3-001 Display
C5	1		33 uF	Capacitor, Tantalum		SBC3-001 Display
C6	0					SBC3-001 Display
C7	0					SBC3-001 Display
C8	1		33 uF	Capacitor, Tantalum		SBC3-001 Display
C9	1		33 uF	Capacitor, Tantalum		SBC3-001 Display
C10	1		1 uF	Capacitor, Tantalum		SBC3-001 Display
C11	1		1 uF	Capacitor, Tantalum		SBC3-001 Display
C12	1		1 uF	Capacitor, Tantalum		SBC3-001 Display
C13	1		1 uF	Capacitor, Tantalum		SBC3-001 Display
C14	0					SBC3-001 Display
C15	0					SBC3-001 Display
C16	0					SBC3-001 Display
C17	0					SBC3-001 Display
C18	0					SBC3-001 Display
C19	0					SBC3-001 Display
C20	0					SBC3-001 Display
C21	0					SBC3-001 Display
C22	0					SBC3-001 Display
C23	0					SBC3-001 Display
C24	0					SBC3-001 Display
C25	0					SBC3-001 Display
C26	0					SBC3-001 Display
C27	0					SBC3-001 Display
C28	0					SBC3-001 Display
CX1	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX2	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX3	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX4	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX5	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX6	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX7	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX8	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX9	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
CX10	1		0.1 uF	Capacitor, mono ceram		SBC3-001 Display
D1	0					SBC3-001 Display
D2	0					SBC3-001 Display
D3	0					SBC3-001 Display
D4	0					SBC3-001 Display
D5	1	1N4148		Diode, fast switching		SBC3-001 Display
D6	1	1N4148		Diode, fast switching		SBC3-001 Display

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
D7	1	1N4148		Diode, fast switching		SBC3-001 Display
D8	1	1N4148		Diode, fast switching		SBC3-001 Display
D9	1	1N4148		Diode, fast switching		SBC3-001 Display
D10	1	1N4148		Diode, fast switching		SBC3-001 Display
D11	1	1N4148		Diode, fast switching		SBC3-001 Display
D12	1	1N4148		Diode, fast switching		SBC3-001 Display
D13	1	1N4148		Diode, fast switching		SBC3-001 Display
D14	1	1N4148		Diode, fast switching		SBC3-001 Display
D15	1	1N4148		Diode, fast switching		SBC3-001 Display
D16	1	1N4148		Diode, fast switching		SBC3-001 Display
D17	1	1N4148		Diode, fast switching		SBC3-001 Display
D18	1	1N4148		Diode, fast switching		SBC3-001 Display
D19	1	1N4148		Diode, fast switching		SBC3-001 Display
D20	1	1N4148		Diode, fast switching		SBC3-001 Display
D21	0					SBC3-001 Display
D22	0					SBC3-001 Display
D23	0					SBC3-001 Display
D24	0					SBC3-001 Display
D25	0					SBC3-001 Display
D26	0					SBC3-001 Display
D27	0					SBC3-001 Display
Q1	1	IRF530		FET, Power, N channel		SBC3-001 Display
U1	1	74HC573		Octal Latch		SBC3-001 Display
U2	1	24C02		EEPROM		SBC3-001 Display
U3	1	80C31BH		8 Bit Microprocessor		SBC3-001 Display
U4	1	MAX813L		Supervisor/Watchdog	Maxim	SBC3-001 Display
U5	1	DMC-16202NY-LY		LCD 2 x 16 Supertwist W/Backlight	Optrex	SBC3-001 Display
U6	1	27C256 W/Program		EPROM		SBC3-001 Display
U7	0					SBC3-001 Display
U8	1	74HC139		Address Latch		SBC3-001 Display
U9	1	74HC00		Quad NAND Gate		SBC3-001 Display
U10	1	LM2940-5.0		Volt Regulator	National Semi	SBC3-001 Display
U11	1	ICL7660		Volt Generator		SBC3-001 Display
U12	1	MAX232		RS232 Driver/PS	Maxim	SBC3-001 Display
U13	0					SBC3-001 Display
U14	0					SBC3-001 Display
U15	0					SBC3-001 Display
U16	1	DS2003		Buffer	National Semi	SBC3-001 Display
U17	0					SBC3-001 Display
U18	1	74HC574		Octal Latch		SBC3-001 Display
U19	1	74HC573		Octal Latch		SBC3-001 Display
U20	1	4N32P		Optocoupler		SBC3-001 Display
U21	0					SBC3-001 Display
U22	0					SBC3-001 Display
U23	0					SBC3-001 Display
U24	0					SBC3-001 Display
U25	0					SBC3-001 Display
U26	1	4013B		Flip Flop		SBC3-001 Display

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
U27	0					SBC3-001 Display
U28	0					SBC3-001 Display
S1	1	EVQ		Switch, Momentary		SBC3-001 Display
K1	1	G5V-1		Relay, SPDT	Omron	SBC3-001 Display
X1	1		7.3728 MHz	Crystal		SBC3-001 Display
X2	0					SBC3-001 Display

DETECTOR VFC BOARD

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
R1	0			Resistor, 1/4W, 5%, CF		6030U-001
R2	0			Resistor, 1/4W, 5%, CF		6030U-001
R3	1		10 Ohm	Resistor, 1/4W, 5%, CF		6030U-001
R4	0			Resistor, 1/4W, 5%, CF		6030U-001
R5	1		10 Ohm	Resistor, 1/4W, 5%, CF		6030U-001
R6	0			Resistor, 1/4W, 5%, CF		6030U-001
R7	1		100K	Trimmer, 20T, Top, 3/8		6030U-001
R8	1		8.06K	Resistor, 1/4W, 1%, MF		6030U-001
R9	0			Resistor, 1/4W, 5%, CF		6030U-001
R10	0			Resistor, 1/4W, 5%, CF		6030U-001
R11	0			Resistor, 1/4W, 5%, CF		6030U-001
R12	0			Resistor, 1/4W, 5%, CF		6030U-001
R13	0			Resistor, 1/4W, 5%, CF		6030U-001
R14	1		8.06K	Resistor, 1/4W, 1%, MF		6030U-001
R15	1		100K	Trimmer, 20T, Top, 3/8		6030U-001
R16				Resistor, 1/4W, 5%, CF		6030U-001
R17	1		10K	Resistor, 1/4W, 5%, CF		6030U-001
R18	1		1.5K	Resistor, 1/4W, 5%, CF		6030U-001
R19	0			Resistor, 1/4W, 5%, CF		6030U-001
R20	0			Resistor, 1/4W, 5%, CF		6030U-001
R21	0			Resistor, 1/4W, 5%, CF		6030U-001
R22	1		2K	Resistor, 1/4W, 5%, CF		6030U-001
R23	1		2K	Resistor, 1/4W, 5%, CF		6030U-001
R24	0			Resistor, 1/4W, 5%, CF		6030U-001
R25	1		10 Ohm	Resistor, 1/4W, 5%, CF		6030U-001
R26	1		100K	Trimmer, 20T, Top, 3/8		6030U-001
R27	1		100K	Resistor, 1/4W, 5%, CF		6030U-001
R28	0			Resistor, 1/4W, 5%, CF		6030U-001
R29	0		100 Ohm	Resistor, 1/4W, 5%, CF		6030U-001
R30	1	MOX300	200M	Resistor, 1/4W, 5%, CF	Victoreen	6030U-001
R31	1		430K	Resistor, 1/4W, 5%, CF		6030U-001
R32	1		2.7K	Resistor, 1/4W, 5%, CF		6030U-001
R33	1		1K	Resistor, 1/4W, 5%, CF		6030U-001
R34	1	MOX300	200M	Resistor, 1/4W, 5%, CF	Victoreen	6030U-001
R35	0			Resistor, 1/4W, 5%, CF		6030U-001
R36	1		499K	Resistor, 1/4W, 1%, MF		6030U-001

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
R37	1		100K	Trimmer, 20T, Top, 3/8		6030U-001
R38	1		51K	Resistor, 1/4W, 5%, CF		6030U-001
R39	1		100K	Resistor, 1/4W, 5%, CF		6030U-001
R40	0			Resistor, 1/4W, 5%, CF		6030U-001
R41	0			Resistor, 1/4W, 5%, CF		6030U-001
R42	0			Resistor, 1/4W, 5%, CF		6030U-001
R43	0			Resistor, 1/4W, 5%, CF		6030U-001
R44	0			Resistor, 1/4W, 5%, CF		6030U-001
R45	0			Resistor, 1/4W, 5%, CF		6030U-001
R46	0			Resistor, 1/4W, 5%, CF		6030U-001
R47	0			Resistor, 1/4W, 5%, CF		6030U-001
R48	1		10M	Resistor, 1/4W, 5%, CF		6030U-001
R49	1		100 Ohm	Resistor, 1/4W, 5%, CF		6030U-001
R50	1		100K	Resistor, 1/4W, 5%, CF		6030U-001
R51	1		100K	Resistor, 1/4W, 5%, CF		6030U-001
C1	0					6030U-001
C2	1		0.1 uF	Capacitor, mono ceram		6030U-001
C3	1		10 uF	Capacitor, Tantalum		6030U-001
C4	1		0.1 uF	Capacitor, mono ceram		6030U-001
C5	1		0.1 uF	Capacitor, mono ceram		6030U-001
C6	1		0.1 uF	Capacitor, mono ceram		6030U-001
C7	1		0.1 uF	Capacitor, polypropelyne		6030U-001
C8	1		0.022 uF	Capacitor, polypropelyne		6030U-001
C9	1		0.01 uF	Capacitor, polypropelyne		6030U-001
C10	1		0.01 uF	Capacitor, mono ceram		6030U-001
C11	0					6030U-001
C12	0					6030U-001
C13	1		0.1 uF	Capacitor, mono ceram		6030U-001
C14	1		0.1 uF	Capacitor, mono ceram		6030U-001
C15	0					6030U-001
C16	1		1000 uF	Capacitor, Electrolytic		6030U-001
C17	1		0.1 uF	Capacitor, mono ceram		6030U-001
C18	1		0.1 uF	Capacitor, mono ceram		6030U-001
C19	0					6030U-001
C20	1		100 uF	Capacitor, Tantalum		6030U-001
C21	1		0.1 uF	Capacitor, mono ceram		6030U-001
C22	1		100 uF	Capacitor, Tantalum		6030U-001
C23	1		0.022 uF 3kV	Capacitor, Disc Ceram		6030U-001
C24	1		100 uF	Capacitor, Tantalum		6030U-001
C25	1		0.022 uF 3kV	Capacitor, Disc Ceram		6030U-001
C26	1		0.022 uF 3kV	Capacitor, Disc Ceram		6030U-001
C27	1		0.022 uF 1.6kV	Capacitor, Polyproplene		6030U-001
C28	1		0.022 uF 3kV	Capacitor, Disc Ceram		6030U-001
C29	0					6030U-001
C30	1		0.001 uF	Capacitor, mono ceram		6030U-001

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
C31	1		100 pF	Capacitor, mono ceram		6030U-001
C32	0					6030U-001
C33	1		0.1 uF	Capacitor, mono ceram		6030U-001
C34	0					6030U-001
C35	0					6030U-001
C36	0					6030U-001
C37	1		0.05 uF 1.6kV	Capacitor, Polypropylene		6030U-001
						6030U-001
CX1	1		0.1 uF	Capacitor, mono ceram		6030U-001
CX2	1		0.1 uF	Capacitor, mono ceram		6030U-001
						6030U-001
D1	0					6030U-001
D2	0					6030U-001
D3	0					6030U-001
D4	0					6030U-001
D5	1	FR107	1 kV	Rectifier, Fast recovery		6030U-001
D6	1	FR107	1 kV	Rectifier, Fast recovery		6030U-001
D7	1	FR107	1 kV	Rectifier, Fast recovery		6030U-001
D8	1	FR107	1 kV	Rectifier, Fast recovery		6030U-001
D9	0					6030U-001
D10	0					6030U-001
D11	0					6030U-001
D12	0					6030U-001
T1	1	6030-T1		Transformer	HPI	6030U-001
U1	1	TLC271CP		Op Amp		6030U-001
U2	1	VFC121		Volt to Freq Converter	Burr Brown	6030U-001
U3	1	75176		RS485 Line Driver/Rcvr		6030U-001
U4	1	4538		Dual Mono		6030U-001
U5	1	LM2940		Volt Reg		6030U-001
U6	1	LM336-2.5		Volt Reference 2.5 V		6030U-001
U7	1	LM336-2.5		Volt Reference 2.5 V		6030U-001
U8	1	ICL7660		Volt Converter		6030U-001
U9	1	TLC271CP		Op Amp		6030U-001
U10	1	ICL7665		Dual Timer		6030U-001
U11	1	TLC27M2		Dual Op Amp		6030U-001
U12	0					6030U-001
U13	0					6030U-001
U14	0					6030U-001
U15	0					6030U-001

PREAMP ELECTROMETER BOARD

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
R101	1	104	5 x 10E10 Ohm	Hi Meg Resistor	Eltec	6030-006
R102	1		10 Ohm	Resistor, 1/4W, 5%, CF		6030-006
R103	1	MOX300	1000M	Resistor, 1/4W, 5%, CF	Victoreen	6030-006
R104	0					6030-006

R105	0					6030-006
R106	1		10 Ohm	Resistor, 1/4W, 5%, CF		6030-006
R107	1		100.0K	Resistor, 1/4W, 1%, MF		6030-006
R108	0					6030-006
R109	0					6030-006
R110	0					6030-006
C101	1		0.1 uF	Capacitor, Mono Ceram		6030-006
C102	1		100 pF	Capacitor, Polystyrene		6030-006
C103	1		10 uF	Capacitor, Tantalum		6030-006
C104	1		1000 pF	Capacitor, Polystyrene		6030-006
C105	1		0.1 uF	Capacitor, Mono Ceram		6030-006
C106	1		0.1 uF	Capacitor, Mono Ceram		6030-006
C107	1		10 uF	Capacitor, Tantalum		6030-006
U101	1	LMC6041		Op Amp	National Semi	6030-006
U102	0	4052B		SS Switch		6030-006
Q101	0					6030-006

DETECTOR CHASSIS

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
M1	1	6030-M1		Chassis Back Cover		
M2	1	6030-M2		Chassis Surround		
M3	1	6030-M3		Chassis Front Cover		
M4	1	6030-M4		Electrometer Box with Cover		
M5	1	1055-6030		Detector		
M6	1	6030-M6		Detector Mounting		
M7	1	6030-M7		Detector Straps		
M8	1	6030-M8		Electrometer Circuit Board Complete		
M9	1	6030-M9		VFC Circuit Board Complete		
M10	1	6030-M10		Cable Electrometer to VFC		
M11	1	6030-M11		Cable VFC to Circular Chassis		
M12	1	6030-M12		Desiccant		
M13	1	6030-M13		Mounting Bracket		
M14	2	6030-M14		Thumbscrews		

READOUT CHASSIS

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
M1	1	6025-M1		Chassis Back Cover		
M2	1	6025-M2		Chassis Bottom Cover		
M3	1	6025-M3		Chassis Front Cover		
M4	1	6025-M4		Chassis Top Cover		
M5	1	6025-M5		Power Supply w/MOV and Jumpers for 110VAC		
M6	1	6025-M6		SBC3 Large Circuit Board		
M7	1	6025-M7		LCD Display		

DESIGNATOR	QUAN	PART NUMBER	VALUE	DESCRIPTION	MFG	DRAWING NO.
M8	1	6025-M8		Push Button Circuit Board		
M9	1	6025-M9		Buzzer		
M10	1	6025-M10		Line Cord		
M11	1	6025-M11		Line Recepticle w/Fuse Holder		
M12	1	6025-M12		Power Switch		
M13	1	6025-M13		Rear Panel I/O Connector		
RX1	1		100 Ohm	Resistor 5W		
D1	1	6025-M14	Red	LED Annunciator w/Circuit Board		
D2	1	6025-M15	Yellow	LED Annunciator w/Circuit Board		
D3	1	6025-M16	Green	LED Annunciator w/Circuit Board		
Dx blank	1	6025-M17		BLANK LED Annunciator wo/Circuit Board		
M14	1	6025-M18		Wiring Harness w/Connectors		
M15	1	6025-M15		Cable LCD to Display PCB		
M16	1	6025-M16		Cable Push Button to Display PCB		
M17	1	6025-M17		Cable LED Annunciators to Display PCB		
M18	1	6025-M18		Cable 6025 to 6030 6' Long		
F1	1	6025-M19		Fuse 1/2 A 5 x 20 mm		

XIII. SCHEMATICS

The Schematics are on the following pages.